

**MOBILE OPERATING TABLE  
ALPHASTAR PRO**





## THE ALPHASTAR PRO OR TABLE – EXCEPTIONAL STABILITY AND UNIVERSAL APPLICATIONS MAQUET – THE GOLD STANDARD



**Always built on practical experience, for even more efficient workflows:** MAQUET is a subsidiary of the Swedish, stock market-listed GETINGE Group. The company is one of the leading global suppliers of therapeutic treatments as well as medical products and OR, catheterisation laboratory and intensive care unit solutions. With its three divisions – Cardiovascular, Critical Care und Surgical Workplaces – MAQUET stands for innovation and technological progress.

**Modular, versatile and extremely stable in use:** The mobile ALPHASTAR PRO OR table. Like all the other products in the ALPHASTAR OR table range, the outstanding features of the ALPHASTAR PRO are its exceptional column stability and the low-level version for operations carried out from the seated position. New: The motorised leg plate mounting point for enhanced flexibility and increased application versatility.

MAQUET – The Gold Standard.

## AS FLEXIBLE AS THE DAILY CLINICAL REQUIREMENTS AS SIMPLE AS THE OR TEAM WANTS

### **Maximum modularity thanks to the use of uniform**

**mounting points:** The use of uniform mounting point geometry means that ALPHASTAR PRO offers the ideal prerequisites for a wide range of positioning options. In the gynaecology and urology disciplines, for example, optimum positioning can be achieved by the use of motor-adjustable leg plates and leg holders.

### **Easy Click for quick changes:**

Quick and simple table-top configuration is made easy thanks to the specially designed mounting point. One single action is all it takes to exchange the motor-adjustable leg and back plates. An audible "click" signals that the module is securely locked into position. No need for time-consuming tightening or turning. Efficient and economic.



*Universally applicable: The mobile OR table with motorised leg plate mounting point*



*Quick and efficient: The modular construction with mounting points for extension modules, adaptable to each surgical discipline*



**Two versions for maximum height adjustability:** As well as the normal version with height adjustment from 685 to 1120 mm (27 to 44"), ALPHASTAR PRO also has a lower height version, from 598 to 950 mm (24 to 37"). So, depending on circumstances and personal preference, operations can be carried out in an equally ergonomic and relaxed way, whether seated or standing.

**Backlit hand control:** Working in a darkened MIS theatre is made much easier thanks to the backlit keypad on the hand control. The icons are self-explanatory and give the user intuitive and secure operation via an IR signal or a cable.



*Choose to work standing up or sitting down: ALPHASTAR PRO is available in the normal or lower height version, depending on your requirements*



*Adjusting height, tilt, slope, leg and back plates using the backlit hand control*

## VERSATILE, ERGONOMIC, COMFORTABLE PROPERTIES THAT COUNT

**Universally applicable:** From the autodrives to the exceptionally stable column with its three-point base, from the electro-hydraulic drive for all the important functions of the OR table to the modular construction of the table top and the mounting point for motorised leg plates, ALPHASTAR PRO can be deployed for all surgical disciplines.

**Intraoperative imaging:** The complete table top needs no crossbars, making it very suitable for intraoperative radioscopes and providing the best conditions for efficient workflows and good operation results.

*Electrical autodrives with smooth start-up and safety brake function, as option with the standard and lower height version*

*Maximum permitted patient weight for the ALPHASTAR PRO: 454 kg / 1000 lbs*

*Modular design for ideal adaptation to a wide variety of surgical tasks*





*Extra thick (80 mm / 3") and tried and tested: Multi-layer SFC pads offer patients protection against pressure injuries and deliver noticeably more comfort*



*Even more safety: The safety catch for the head plate mounting point (passive lock)*



*Active unlocking when removing the head plate*

*Option: Joint adapter (1131.82A0) for more accessories*

*Ergonomic working position, whether sitting or standing, thanks to flexible height adjustability*

*Enhanced modularity: Mounting point for extension modules*



## GYNAECOLOGICAL AND UROLOGICAL POSITIONING MOTORISED LEG PLATES FOR ADDED PRECISION

**Optimum positioning best results:** ALPHASTAR PRO offers the best conditions for working comfort with height adjustments to meet all expectations, motor-adjustable leg plates and its narrow column foot for ease of access. The seat plate

extension with leg holders facilitates the complete radioscopy of the urogenital tract. The leg plate mounting point can be moved upwards through a maximum of 80° and downwards through 90°.



*Lithotomy position using Goepel knee crutches*



*Lithotomy position using Goepel knee crutches*



*Lithotomy position with seat plate extension and leg holders with single-handed operation*



*Lateral position for kidney operations*

\* radiolucent area





## ORTHOPAEDIC AND TRAUMATOLOGY POSITIONING VERSATILE EXPANSION OPTIONS AND 360° RADIOTRANSLUCENT

**The best conditions for every type of operation:** Uniform mounting points and simple "Easy Click" installation make quick work of equipping ALPHASTAR PRO with MAQUET special orthopaedic and traumatology accessories.

Intraoperative controls via imaging procedures are made possible, thanks to the 360° radiotranslucent modules, giving optimum results and reducing the need for post-operative measures.



*Beach chair position with TRIMANO 3D support arm*



*Lateral position for hip operations*



*Prone position for spinal operations*



*Prone position with positioning device for spinal surgery, giving maximum stability and optimum a-p and 360° radioscopy of the spine*

## POSITIONING EXAMPLES FURTHER DISCIPLINES



*"Reverse Mode": - 360° fluoroscopy possible with image converter (C-arm) for surgery to lower extremities or pelvis.*



*Gentle head support, maintaining the correct spinal position, using a horseshoe-shaped head rest with motor-driven head rest adjustment*




*"Normal Mode": Unrestricted access for cardiovascular surgery with the image converter (C-arm), from the head right down to the lower leg*



*Position for laparoscopic adiposity surgery*

## ALPHASTAR PRO ACCESSORY OVERVIEW





*Seated position  
for neurosurgery*



## TECHNICAL SPECIFICATIONS AND CONSTRUCTION FEATURES OF ALPHASTAR PRO

### Adjustment options using the cabled control unit, IR control unit or override panel

Height without padding

1132.21

685-1120 mm / 27-44"

1132.22

598-950 mm / 24-37"

Trendelenburg

+30°/-30°

Lateral tilt, left/right

18°

Normal back plate

+70°/-40°

Leg plate mounting point

+80°/-90°

0-position (resetting of adjustment

options such as slope, lateral tilt, back

and leg plates)

Base foot locking (Lock/Unlock)

Patient orientation (normal/reverse)

Base foot autodrive

(Equipment option, in conjunction with corded hand control)

### Technical information

Max. patient weight

1132.21/22

454 kg / 1000 lbs

Operating table weight

300 kg / 661 lbs

CE Complies with Medical Device Directive 93/42/EEC

### Dimensions

Length without accessories

860 mm / 34"

Width at the back plate










525 mm / 21"

Width at the seat plate

540 mm / 21"

Width across accessory rails

590 mm / 23"



## GENERAL CONSTRUCTION FEATURES ALPHASTAR PRO

- Battery and mains power operation (see electrical specifications)
- Control of the electro-hydraulic manoeuvrability and adjustment functions using a corded hand control or IR remote control
- Stable base construction with four double-swivel castors for easy movement and manoeuvring (locking via the control unit)
- Optional additional foot-operated switch with all models
- Cover for base and override control panel made of glass fibre-reinforced composite plastic, resistant to impact, breakage and disinfectant agents, with neutral tone colouring and a scratch-resistant coating
- CrNi steel column casing
- CrNi steel seat and back plate supporting arms
- CrNi steel accessory rails

## FEATURES OF THE OPERATING TABLE TOP:

- Operating table top: Subdivided into six sections with head plate (optional), upper back plate (optional), lower back plate, seat plate, leg plates (optional)
- Seat and back plate supporting arms, leg clamps, joint covers and accessory rails made of CrNi steel
- SFC (Soft Foam Core) padding with multi-layer construction, for exceptional positioning comfort and maximum patient safety, with 80 mm (3") pad thicknesses
- Table top without crossbars for intraoperative radioscopy
- Setting of patient orientation using the control units

## ELECTRICAL SPECIFICATIONS:

- Maintenance-free battery with capacity for about one week of OR use
- Battery charge level: Electronic monitoring with optical and acoustic display
- Power cable for mains power operation and battery charging
- AC mains power nominal voltages  
100/110-115/127/200/220/230-240 V (switchable)  
Protection class II, Type B, IEC 60601-1  
Rated frequency 50/60 Hz

## VERSIONS AND TABLE TOP PADDING:

1132.21B0	<b>ALPHASTAR PRO Standard version</b> , max. patient weight 454 kg (1000 lbs). Mobile operating table for general surgery in modular design, electro-hydraulic, without control units, with 80 mm (3") SFC padding, radiotranslucent and electrically conductive, 10 x 25 mm (0.3 x 0.9") format accessory rails
1132.21B2	see above (1132.21B0), however, with electric autodrive
1132.22B0	<b>ALPHASTAR PRO low height version</b> , max. patient weight 454 kg (1000 lbs)
1132.22B2	see above (1132.22B0), however, with electric autodrive

## RECOMMENDED BASIC ACCESSORIES:

1132.90J0	Corded hand control
1130.64C0	Tilttable head plate, with SFC padding
1131.31BC	Standard back plate for general surgery, with SFC padding
1133.53BC	Pair of leg plates, with dual-joint abduction, with SFC padding

## OPTIONAL ACCESSORIES:

### Operation via:

1132.91J0	IR control unit
1009.70A0	Charger for IR control unit
1009.81E0	Foot switch for the following functions: height, slope, back
1009.81E1	Foot switch for the following functions: height, slope, tilt
1009.81E2	Foot switch for the following functions: height, slope, leg plates

### Head plates:

1130.81A0	Head plate adapter
1131.82A0	Joint adapter, for 1133 or 1150 modules
1130.67B0	Head plate with gas strut assisted adjustment, with SFC padding
1130.53B0	Head plate with double joint, with SFC padding

### Back/seat plates:

1132.34A0	Back plate for shoulder operations, with SFC padding
1131.55BC	Seat plate extension, with SFC padding

### Leg plates:

1133.58BC	Single-section leg plate, with SFC padding
1132.65A0	Leg positioning device, max. 40 kg (88 lbs) load (weight of a leg from a 135 kg (297 lbs) patient), without accessory rails

**For local contact:**

Please visit our Website  
[www.maquet.com](http://www.maquet.com)

**MAQUET**  
GETINGE GROUP

MAQUET GmbH  
Kehler Str. 31  
76437 Rastatt, Germany  
Phone: +49 7222 932-0  
Fax: +49 7222 932-571  
[info.sales@maquet.de](mailto:info.sales@maquet.de)  
[www.maquet.com](http://www.maquet.com)

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.