# **Technical Specifications**

#### **Temperature Monitor**

Display 30° to 100°C in all lesion/pulsed modes

Temperature Probes Wide range of NeuroTherm probes available

Impedance Monitor **Display** Reads biological impedance in all different modes throughout the procedure

Range 50-2000 Ohms (in 1 Ohm steps)

Self test Internal 500 Ohm test resistor

Sound Audible impedance monitoring available

#### **Stimulation**

Waveform Balanced biphasic waveform with negative leading edge

Pulse width 0.1, 0.2, 0.5, 1.0 mSec

Pulse rate Sensory 10, 20, 50, 75, 100, 150, 180, 200 Hz Motor 2, 5 Hz

## Amplitude Constant Voltage

0-5.0 V 0-3.0 V 0-0.5 V 0-10 mAmp 0-6 mAmp Constant Current 0-1 mAmp

Motor and sensory ranges are independent of each other

**Continuous RF Lesion Mode** RF power 0-30 Watts into 200 Ohms

Frequency 480 kHz sine wave

Display Voltage 0-99 Volts (RMS) Current 0-999 mAmps (RMS) Power 0-30 Watts Impedance 50-2000 Ohms Dynamic graphic display of temperature 30-100°C

Time Selectable 0:30 to 10:00 Timer automatically starts when temperature is within 5°C of  $$-10^{-1}$$ 

Auto mode

By pressing the Auto start button the temperature will ramp up 8°C per second until set temperature is reached. The time will start when the temperature is within 5°C

Maximum temperature set Selectable 50-90°C in 5° or 1°C steps. Automatically adjusts RF power to not exceed temperature.

#### IDET

P1, P2 and P3 are three pre-programmed IDET Profiles. In custom step mode there is the possibility to set up different IDET profiles by setting Start temp, Step time, Step rise, Final temp. Final dwell time

#### **Pulsed RF Mode**

Output User definable pulse bursts of RF power

## Pulse burst width 5, 10, 20, 50 mSec

Pulse burst frequency 1, 2, 5, 10 Hz (Only 1, 2 Hz in multiple Probe Mode)

Amplitude range Constant Voltage 30 - 70 Volts (RMS) 50 - 350 mAmps (RMS) Constant Curren

Maximum temperature set Selectable in 42 – 90°C range

Selectable 0:30 – 20:00 minutes

#### Pulsed Dose In pulsed dose mode only pulses of the full pulse width and full amplitude will be delivered to the patient. The number of pulses to be delivered is chosen in place of time. A counter displays the delivered dose. 120 to 2400 Pulses

#### **Multiple Probes**

The NeuroTherm NT1100 can be used with 1, 2 or 3 probes simultaneously in RF lesion, Pulsed RF and Pulsed Dose modes. The "Intellipower" algorithm ensures a smooth and safe. synchronous temperature rise of each electrode

**Safety Features** Sterile probe test Checks probe intraoperatively

#### Output lockout

Hardware and software lockout if voltage / current control is not initially set to zero in all stimulation and RF modes. Maximum temperature Lesion temperature is limited to 90°C Additional hardware lockout should temperature exceed 95°C

Thermocouple Lockout when a faulty thermocouple electrode is connected to one of the sockets

#### Operation

All settings are entered through the touch screen. When actually performing the procedure the touch screen is only used as display. All active controls are located on the control panel.

#### Data Input User defined settings Up to 12 different doctors can store their own unique settings

in the internal memory of the generator. Patient details Patient details can be entered into internal memory and used in reports. After entering the patient details, the patients can be selected in a drop-down menu, when starting the procedure

#### **Data Acquisition**

Data transfer All procedure details can be exported out of the NT1100 in CSV file format. These files can be converted into for example Microsoft ExcelĐ files.

#### Instant procedure reports

Reports of all procedures performed are available on the NT1100. These reports can be viewed on the screen, printed on the Bluetooth printer or exported to the USB stick in TXT file format

**Procedure site labelling** Up to eight custom "site-labels" can be entered into the NT1100, to indicate the different treatment sites. Site labels can be assigned to a treatment by a simple touch on the screen, before or during the procedure

NeuroTherm<sup>®</sup> 30 Upton Drive, Suite 2 Wilmington, MA 01887 USA

Tel: +1 888 655 3500 Fax:+1 978 658 2378

## EC REP

**Peripheral Accessories** 

connection, to print procedure reports.

computer

Additional Screen

Defaults

80°C 1.00 min

Pulsed RF

Pulse Dose

Standards EN60601-1:1997

EC60601-1-21983

IEC60601-2-2:1998

Warranty

corporation, USA

Two years

Sensory stimulation 50 Hz, 1.0 mSec, 03 V

Motor stimulation 2 Hz, 1.0 mSec, 03 V

Continuous RF lesion mode

20 mSec, 2 Hz, 42°C, 45 V, 2:00 min

20 mSec, 2 Hz, 42°C, 45 V, 240 pulses

IEC60601-2-10: With Canadian deviations

C22.2 No. 6011 and IEC60601-1

Protection for electrical shock Class II type BF

With respect to electrical shock, fire and mechanical hazards only in accordance with UL60601-1, IEC60601-1, CAN/CSA

Any equipment connected to rear sockets must comply with IEC60950 and IEC60601-1

Bluetooth™ is a registered trademark of Bluetooth SIG, Inc.,

USA Microsoft Excel™ is a registered trademark of Microsoft

USB connection A USB memory stick (delivered with the NT1100) can be

the NT1100 with the NeuroTherm video unit (NT-VD)

used to export the procedure details from the NT1100 to a

An external screen can be connected to the NT1100 giving a repeat view of the display screen. The output is isolated from

**Printer** A printer can be connected to the NT1100 using a Bluetooth

**NeuroTherm**<sup>®</sup> 429 Brighton Road, Croydon Surrey, CR2 6EU United Kingdom Tel: +44 20 8660 4374 Fax:+44 20 8660 9417

**NeuroTherm®** Modemstraat 20 1033 RW Amsterdam The Netherlands Tel: +31 204 702 702 Fax:+31 204 702 701

# **NeuroTherm**<sup>®</sup>



www.neurotherm.com

Distributed by:

# ISO 13485 Certified

**Always** a step ahead...



The first multi-Electrode **RF** Generator now has even more to offer.

NT 1100, the newest and most advanced Radiofrequency Generator, is the only machine compatible with the latest technologies designed by NeuroTherm. This software based generator sets the standard for innovation, efficiency, and safety.

# Customize your NT 1100 Radiofrequency Generator.

Program your personal preferences for every mode or procedure, add patient data at the start of your day or prior to their procedure, and label up to 8 anatomical sites.

NT 1100 comes with integrated data acquisition software so you can access your data when you want it.

Data Management right at your fingertips.



Setup Set	election		Frequency Structures Frequency SC Hz Width 1.0 ml Volte 0-3 V
Pulse Dose Set Temp	42 🔻	°C	Votor Standation Frequency 2 Hz Width 1.9 m3 Vots 0-3 V
Pulse Count	and the second second	Ū.	Themail Lesion Set Temp 90 °C Time 2 Mins Step Lesion P1
Rate Width	2 💌 20 💌	PPS mS	Fulsed Set Temp 43 °C Time: 120 Se Rate 2 PPS Width 20 mS
Set Volt	45 💌	v	Putsed Dese Set Temp 42 °C Putse Count 240 Rate 2 PPS With 20 mS Set Yolk 45 V
nerne <sup>ver</sup> ITT 1000 S. R. Senith	23 March 2015	Next	Caution: Machine Stat

Procedure	Set	tting	IS		Frequency 5 Frequency Width Runge	
Electrode Selection	1 +	2	3	Disp +	Motor Star Frequency Width Range	
Alt Single Electrode RF		+		+	Thermal L	-
Two Electrodes RF		+	+	+	Set long	-
Three Electrodes RF	+	+	+	+	Time Step Lesin	1:0
Dual Electrode		+	+	+	Pulsed RF	
Cordotomy Electrode		+		+	Set long	
Bipolar Electrode		+			Time	2:0
Two Bipolar Electrodes		+	+		with	20
No Thermocouple Elect	+				Set	45
					Pulsed Dos	
					Set Temp Court	42
					Rate	29
		100			widh	28
			Acce	pt	Set	45

#### Pulse Dose\*

Finish

The alternative that makes sense. For the first time ever, you can be sure of the voltage delivered with every pulse of radiofrequency. Using Pulse Dose will not compromise amplitude, pulse width, or pulse rate. Simply set your preferred voltage, press auto start, and count down the pulses.

Doctor Settings

2 3 4 5 6 7 8 9 0 Q W E R T Y U I O P

Date of Birth Referring Doct dditional Detail

Delete All

Site 4

Site 5

Site 6

Edit

Accept

20/Sep/2007 20/Sep/2007 20/Sep/2007 21/Sep/2007 12/Nov/2007 12/Nov/2007

New

Delete Selected Copy To USB Stick Print Finish

ZXCVBNM

Site Labels

4 5 6 7 8 9 0 R T Y U I O F

HJKL

File Control View Print Files Copy CSV files

Copy BMP files

Delete CSV files Delete BMP files

Print session data More ..

Accept

Site 1 Description L3-L4 LEFT

ZXCVBNM

#### Dual Electrode

Only the NT 1100 has a dual electrode capability for Thermal Lesion and Pulsed Radiofrequency modes where RF energy is communicated between the two active tips. Optimal placement of the two electrodes can create a larger lesion size or electric field compared to using only one electrode. A safe and economical clinical option.

### Multiple Electrodes

The first to offer Multiple Electrodes in simultaneous use. Perform RF procedures more efficiently by treating up to 3 sites simultaneously in either Lesion, Pulsed, or Pulse Dose modes. Each electrode is controlled independently to maximize safety during your procedure.

# Simplicity **\***

Simplicity III is a new probe designed with three distinct active areas. It is flexible enough to navigate challenging anatomical areas and the user can control the shape of lesion or electric field created in the nerve tissue. Harness the power of RF with Simplicty III.

# Diskit II

Diskit II offers a safe, simple, and effective alternative to disc treatments. The kit comprises of two easy to place 20 gauge introducers with matching electrodes. The unique algorithm in the NT 1100 directs the RF energy between the two active tips without the use of a dispersive plate. Diskit II can be used either in Thermal lesion or Pulsed RF to treat discogenic pain. The result is a larger lesion within the annulus or a wider electric field in Pulsed mode.

www.neurotherm.com