

Phlebotomy Chair with Hi-Lo and Power Back

Model 8678


The UMF Medical 8678 provides Power Hi-Lo & Power Backrest movement with a two-function hand control with OneTouch Patient Positioning® Reclines to 30 degrees from horizontal position with self-leveling armrests. Meets ADA height requirement and has a 375lbs weight capacity.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



Key features at a glance

- › Ultra Comfort Seamless Top with antimicrobial guard
- › 19" (48.3cm) ADA compliant low access height
- › Power Hi-Lo
- › Power backrest
- › Two function hand control with OneTouch Patient Positioning™
- › Self-positioning armrests
- › Self-positioning leg section
- › Reclines to 30° from horizontal position
- › Electro-mechanical actuator system
- › 375lb (170.1kg) patient weight capacity
- › 3-year limited warranty

Specifications

Seat Height: Min. 19" (48.3cm) at seat and Max. 36" (91.4cm) at seat
 Overall Width: 31.5" (80cm)
 Seat Width: 18" (45.7cm)
 Length (fully extended): 67" (170.2cm)
 Arm Length: 23" (58.4cm)
 Electrical: 120 VAC, 60 HZ or 220/240 VAC, 60HZ
 Backrest Width: 17" Max / 12" Min. (43.18cm Max/ 30.48cm Min.)
 Net Weight: 216lb (98kg)
 Shipping Weight: 275lb (124.8kg)
 Shipping Dimensions: L44xW34xH64 (111.8cm x 86.4cm x 162.6cm)

Accessories

032 – Optional backrest with articulating headrest (factory installed)
 033 – Fold-away table 13.5" W x 10" L x .75" D (34.3cm W x 25.4cm L x 1.9cm H)
 143 – Easy-Glide locking casters (set of 4/factory installed)
 422 – Two-function foot control
 856 – 220/240 Volt (factory installed)
 8678F – Fire retardant foam